



Human IgG protein (Fab'2) [FITC] (DAGA-3143)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Immunoglobulin G protein (Fab'2) (FITC)
Purity	Human IgG protein (Fab'2) (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by papain digestion and dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.01M Na ₃ PO ₄ , pH 7.4, with 0.14M NaCl, 10 mg/ml BSA, and 0.01% thimerisol. Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.