



Norsertraline [BSA] (DAGA-315B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|--|
| Antigen Description | Norsertraline is an antidepressant of the selective serotonin reuptake inhibitor (SSRI) class. |
| Conjugate | BSA |
| Applications | LFIA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|----------|---------------|
| Keywords | Norsertraline |
|----------|---------------|