



Human IgM protein (Fab mu) [FITC] (DAGA-3173)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Human Immunoglobulin Mprotein (Fab mu) (FITC)
Purity	Human IgM protein (Fab mu) (FITC) was purified by delipidation, salt fractionation, ion exchange chromatography and dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA and 0.01% NaN3. Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.