



Guinea Pig IgG (Fab) [HRP] (DAGA-3177)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Guinea Pig Immunoglobulin G (Fab) (HRP)
Purity	Guinea Pig IgG (Fab) (HRP) was purified by delipidation, salt fractionation, ion exchange chromatography followed by dialysis.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.02M K ₃ PO ₄ , pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% gentamicin sulfate. DO NOT add NaN ₃ . Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.