



Turkey IgG protein [HRP] (DAGA-3185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Turkey Immunoglobulin G protein (HRP)
Purity	Turkey IgG protein (HRP) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.02M K ₃ PO ₄ , pH 7.2, with 0.15M NaCl, 10 mg/ml BSA and 0.01% gentamicin sulfate. DO NOT add NaN ₃ . Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.