



Rat IgG (Fab'2) [HRP] (DAGA-3193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Immunoglobulin G (Fab'2) (HRP)
Purity	Rat IgG (Fab'2) (HRP) was purified by delipidation, salt fractionation and ion exchange chromatography.
Conjugate	HRP
Applications	User optimized
Format	Lyophilized from 0.01M Na ₃ PO ₄ , pH 7.4, with 0.14M NaCl, 10 mg/ml BSA and 0.01% thimerisol. Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.