



## Rabbit IgG (Fab) [FITC] (DAGA-3197)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                         |                                                                                                                                                                     |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Rabbit Immunoglobulin G (Fab) (FITC)                                                                                                                                |
| <b>Purity</b>           | Rabbit IgG (Fab) (FITC) was purified by affinity chromatography.                                                                                                    |
| <b>Conjugate</b>        | FITC                                                                                                                                                                |
| <b>Applications</b>     | IF                                                                                                                                                                  |
| <b>Format</b>           | Lyophilized from 0.02M K <sub>3</sub> PO <sub>4</sub> , pH 7.2, with 0.15 M NaCl, 10 mg/ml BSA, and 0.01% NaN <sub>3</sub> . Immunoglobulin and protease free.      |
| <b>Concentration</b>    | Batch dependent - please inquire should you have specific requirements.                                                                                             |
| <b>Size</b>             | 1 mg                                                                                                                                                                |
| <b>Storage</b>          | Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles. |