



Dog IgG (Fab'2) [FITC] (DAGA-3214)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Dog Immunoglobulin G (Fab'2) (FITC)
Purity	Dog IgG (Fab'2) (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.01M Na ₃ PO ₄ , pH 7.4, with 0.14M NaCl, 10 mg/ml BSA, and 0.01% thimerisol. Immunoglobulin and Protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.