



Dog IgG (Fab) [FITC] (DAGA-3225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Dog Immunoglobulin G (Fab) (FITC)
Purity	Dog IgG (Fab) (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.02 M K ₃ PO ₄ , pH 7.2, with 0.15M NaCl, 10 mg/ml BSA, and 0.01% NaN ₃ . Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.