



Phencyclidine [BSA] (DAGA-322B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Phencyclidine modulates secretogranin II[sg II] expression in prefrontal cortex tissue in the absence of afferent inputs. The nature of these changes is dependent upon the duration of exposure to and/or withdrawal from phencyclidine. It is used as an anesthetic by veterinarians and is illicitly taken for its hallucinogen effects.
----------------------------	---

Conjugate	BSA
------------------	-----

Applications	LFIA
---------------------	------

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	0.01M pH7.4 PBS
---------------	-----------------

Preservative	None
---------------------	------

Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---

BACKGROUND

Keywords	1 (1 phenylcyclohexyl)piperidine; Phencyclidine; Phencyclidine
-----------------	--