



## Phencyclidine [KLH] (DAGA-322K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Phencyclidine modulates secretogranin II[sg II] expression in prefrontal cortex tissue in the absence of afferent inputs. The nature of these changes is dependent upon the duration of exposure to and/or withdrawal from phencyclidine. It is used as an anesthetic by veterinarians and is illicitly taken for its hallucinogen effects.
----------------------------	---

<b>Conjugate</b>	KLH
------------------	-----

<b>Applications</b>	Immunogen
---------------------	-----------

<b>Format</b>	Liquid
---------------	--------

<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
----------------------	---

<b>Size</b>	1 mg
-------------	------

<b>Buffer</b>	0.01M pH7.4 PBS
---------------	-----------------

<b>Preservative</b>	None
---------------------	------

<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---

### BACKGROUND

<b>Keywords</b>	1 (1 phenylcyclohexyl)piperidine; Phency; Phencyclidine
-----------------	---