



Phenobarbital [BSA] (DAGA-324B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Phenobarbital, a barbiturate, is used as a sleep aid and in the treatment of certain types of epilepsy, including generalized or grand mal seizures and partial seizures.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Luminal; Phenobarbitone;
----------	--------------------------