



Prednisolone [KLH] (DAGA-326K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Prednisolone is used to achieve prompt suppression of inflammation in many inflammatory and allergic conditions. Examples of inflammatory conditions include rheumatoid arthritis, systemic lupus, acute gouty arthritis, psoriatic arthritis, ulcerative colitis, and Crohns disease. Severe allergic conditions that fail conventional treatment may also respond to prednisolone.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Prednisolone; Pediapred Oral Liquid; Medrol;
-----------------	--