



# Human Albumin protein [HRP] (DAGA-327)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Albumin refers generally to any protein that is water soluble, which is moderately soluble in concentrated salt solutions, and experiences heat denaturation. They are commonly found in blood plasma, and are unique to other plasma proteins in that they are not glycosylated. Substances containing albumin, such as egg white, are called albuminoids.
<b>Purity</b>	Human Albumin protein (HRP) was purified by ion exchange chromatography followed by dialysis.
<b>Conjugate</b>	HRP
<b>Applications</b>	N/A
<b>Format</b>	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml PEG any 0.01% gentamicin sulfate. DO NOT add NaN3
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Preservative</b>	None
<b>Storage</b>	Store at 4 °C until reconstitution. Following reconstitution aliquot and freeze at -20 °C for long term storage. Avoid repeated freeze/thaw cycles.

## BACKGROUND

<b>Keywords</b>	ALB; albumin; FDAH; ANALBA; PRO0883; PRO0903; PRO1341; serum albumin; albumin (32 AA); albumin (AA 34); growth-inhibiting protein 20; cell growth inhibiting protein 42; anti-Serum Albumin
-----------------	---