



Propoxyphene [KLH] (DAGA-327K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dextropropoxyphene, manufactured by Eli Lilly and Company, is an analgesic in the opioid category. It is intended to treat mild pain and has, in addition, anti-tussive and local anesthetic effects. It has been taken off the market in Europe and the US due to concerns of fatal overdoses and arrhythmias. An estimated 10 million patients have used these products.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	(+)-PROPOXYPHENE; PROPOXYPHENE; (+)-1; 2-diphenyl-2-propionoxy-3-methyl-4-dimethylaminobutane; (+)-4-(Dimethylamino)-3-methyl-1; 2-diphenyl-2-butanol propionate;
-----------------	---