



Human Albumin protein [FITC] (DAGA-328)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Albumin refers generally to any protein that is water soluble, which is moderately soluble in concentrated salt solutions, and experiences heat denaturation. They are commonly found in blood plasma, and are unique to other plasma proteins in that they are not glycosylated. Substances containing albumin, such as egg white, are called albuminoids.
Purity	Human Albumin protein (FITC) was purified by ion exchange chromatography followed by dialysis
Conjugate	FITC
Applications	ELISA
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA and 0.01% NaN3. Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Preservative	None
Storage	Store at 4 °C until reconstitution. Following reconstitution aliquot and freeze at -20 °C for long term storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords	ALB; albumin; FDAH; ANALBA; PRO0883; PRO0903; PRO1341; serum albumin; albumin (32 AA); albumin (AA 34); growth-inhibiting protein 20; cell growth inhibiting protein 42; anti-Serum Albumin
-----------------	---