



# Recombinant Acidianus filamentous virus 2 Putative transmembrane protein ORF289 [His] (DAGA-3286)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Antigen Description</b> | Recombinant 6xHis tagged Acidianus filamentous virus 2 Putative transmembrane protein ORF289 (UniProtKB No. Q573C7) was expressed.                                                                                                                                                 |
| <b>Purity</b>              | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                    |
| <b>Conjugate</b>           | His                                                                                                                                                                                                                                                                                |
| <b>Format</b>              | lyophilized                                                                                                                                                                                                                                                                        |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.                                                                                                                                                                                                            |
| <b>Size</b>                | 1 mg                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>              | Tris-based buffer, 50% glycerol                                                                                                                                                                                                                                                    |
| <b>Storage</b>             | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |