



# Recombinant Acidianus two-tailed virus Uncharacterized protein ORF58 [His] (DAGA-3289)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant 6xHis tagged Acidianus two-tailed virus Uncharacterized protein ORF58 (UniProtKB No. Q3V4W8) was expressed.
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	lyophilized
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.