



Recombinant *Acidilobus saccharovorans* tRNA(Ile2) 2-agmatinylcytidine synthetase TiaS [His] (strain DSM 16705), *E. coli* (DAGA-3291)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged <i>Acidilobus saccharovorans</i> tRNA(Ile2) 2-agmatinylcytidine synthetase TiaS (UniProtKB No. D9Q2C4) was expressed.
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.