



Rotenone [KLH] (DAGA-329K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Rotenone is a colorless-to-red, odorless solid. In solution it is used as a broad spectrum insecticide that works by inhibiting the transfer of electrons from Fe-S centers in Complex I to ubiquinone (see electron transfer chain).
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Canex;
-----------------	--------