



# Recombinant *Acidiphilium cryptum* Reaction center protein L chain pufL [His] (DAGA-3300)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |   |
|----------------------------|---|
| <b>Antigen Description</b> | Recombinant 6xHis tagged Reaction center protein L chain pufL (UniProtKB No. O66137) was expressed.   |
| <b>Purity</b>              | >85% (SDS-PAGE)   |
| <b>Conjugate</b>           | His   |
| <b>Format</b>              | lyophilized   |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.   |
| <b>Size</b>                | 1 mg  |
| <b>Buffer</b>              | Tris-based buffer, 50% glycerol   |
| <b>Storage</b>             | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |