



Recombinant *Acidiphilium cryptum* atpB [His] (strain JF-5) (DAGA-3304)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant 6xHis tagged <i>Acidiphilium cryptum</i> atpB) (UniProtKB No. A5FVJ0) was expressed.
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>