



Salvinorin [KLH] (DAGA-331K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Salvinorin A is the main active psychotropic molecule in <i>Salvia divinorum</i> , a Mexican plant which has a long history of use as an entheogen by indigenous Mazatec shamans. Salvinorin A is a hallucinogenic compound with psychedelic/dissociative effects.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Salvinorin; Salvinorin A;