



Salvinorin A [BSA] (DAGA-332B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Salvinorin A is a hallucinogenic compound with psychedelic/dissociative effects. Salvinorin A is one of the most potent naturally occurring psychoactive drugs known to date, with an effective dose in humans in the 200–1,000 µg range when smoked. Salvinorin A can produce psychoactive experiences in humans with a typical duration of action being several minutes to an hour or so, depending on the method of ingestion.

Conjugate BSA

Applications LFIA

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 0.01M pH7.4 PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Salvinorin A; Salvinorin;