



Salvinorin A [KLH] (DAGA-332K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | Salvinorin A is a hallucinogenic compound with psychedelic/dissociative effects. Salvinorin A is one of the most potent naturally occurring psychoactive drugs known to date, with an effective dose in humans in the 200–1,000 µg range when smoked. Salvinorin A can produce psychoactive experiences in humans with a typical duration of action being several minutes to an hour or so, depending on the method of ingestion. |
| Conjugate | KLH |
| Applications | Immunogen |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

Keywords Salvinorin A; Salvinorin;