



## Secobarbital [HRP] (DAGA-333H)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	A short acting barbiturate. Secobarbital sodium is a white odorless slightly bitter powder (trade name Seconal) used as a sodium salt for sedation and to treat convulsions.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	Secobarbital Sodium; Barbosec; component of Hyptran; Evronal; Hypotrol; Imesonol; Immenoctal; Immenox; Meballymal; Quinalbarbital; quinobarbitone; Secobarbitone; Seconal; somosal; Trisomnin; S
----------	--