



Temazepam [KLH] (DAGA-334K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Temazepam is approved for the short-term treatment of insomnia. In addition, temazepam has anxiolytic (anti-anxiety), anticonvulsant, and skeletal muscle relaxant properties.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Temazepam; Restoril; Normison;
----------	--------------------------------