



Tramadol [KLH] (DAGA-338K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tramadol,(1R,2R)-2-(dimethylaminomethyl)-1-(3-methoxyphenyl)cyclohexan-1-ol, is anarcotic analgesic used for severe pain. It may be habituating.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Tramadol; (1R; 2R) 2 (dimethylaminomethyl) 1 (3 methoxyphenyl) cyclohexan 1 ol; Ralivia; Ultram; CRISPIN; cis-Tramadol
----------	--