



## Tricyclic depressants [KLH] (DAGA-342K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Tricyclic antidepressants (TCAs) are heterocyclic chemical compounds used primarily as antidepressants. The TCAs were first discovered in the early 1950s and were subsequently introduced later in the decade; they are named after their chemical structure, which contains three rings of atoms. The tetracyclic antidepressants (TeCAs), which contain four rings of atoms, are a closely related group of antidepressant compounds.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	TCAs; Tricyclic
-----------------	-----------------