



Zaleplon [KLH] (DAGA-344K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Zaleplon is a sedative/hypnotic, mainly used for insomnia. It is a nonbenzodiazepine hypnotic from the pyrazolopyrimidine class. In terms of adverse effects zaleplon appears to offer little improvement compared to both benzodiazepines and other non-benzod
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Zaleplon; Sonata; Starnoc; N-(3-(3-cyanopyrazolo[1; 5-a] pyrimidin-7-yl)phenyl)-N-ethylacetamide; AKOS 92107; CL-284846; Sonata; N-[3-(3-cyanopyrazolo[1; 5-a]py; Zaleplon and intermediates; CL-284846; N-[3-(3-Cyanopyrazolo-[1; 5-a]pyrimidin-7-yl)phenyl)-N-et
-----------------	--