



# Recombinant Protein *Aedes aegypti* mt:ND5, *E.coli* (DAGA-3465)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The following tags are available.

- N-terminal His-tagged and C-terminal Myc-tagged
- N-terminal His-SUMO-tagged and C-terminal Myc-tagged
- N-terminal His-B2M-tagged and C-terminal Myc-tagged
- N-terminal His-GST-tagged and C-terminal Myc-tagged
- N-terminal His-tagged
- N-terminal His-SUMO-tagged
- N-terminal His-B2M-tagged
- N-terminal His-GST-tagged
- N-terminal GST-tagged
- Tag-Free

<b>Antigen Description</b>	Recombinant 6xHis tagged <i>Aedes aegypti</i> mt:ND5 (UniProtKB No. B0FWD3) was expressed.
<b>Species</b>	<i>Aedes aegypti</i>
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.