



Recombinant Protein Aedes albopictus (Asian tiger mosquito, tegomyia albopicta) CECC1, E.coli (DAGA-3469)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description	Recombinant 6xHis tagged Aedes albopictus (Asian tiger mosquito, tegomyia albopicta) CECC1 (UniProtKB No. Q9Y0X9) was expressed.
----------------------------	--

Species	Aedes albopictus
----------------	------------------

Purity	>85% (SDS-PAGE)
---------------	-----------------

Conjugate	His
------------------	-----

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
Buffer	Tris-based buffer,50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.