



Recombinant Protein Aedes albopictus (Asian tiger mosquito, tegomyia albopicta) CECC1, Baculovirus (DAGA-3470)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged

N-terminal His-tagged

Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description

Recombinant 6xHis tagged Aedes albopictus (Asian tiger mosquito, tegomyia albopicta) CECC1 (UniProtKB No. Q9Y0X9) was expressed.

Species

Aedes albopictus

Purity

>85% (SDS-PAGE)

Conjugate

His

Format

Liquid

Concentration

Batch dependent - please inquire should you have specific requirements.

Size

1 mg

Buffer

Tris-based buffer, 50% glycerol

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.