



Recombinant Protein *Aeromonas hydrophila* subsp. *Hydrophila* nqrD, Mammalian cell (DAGA-3501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged <i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i> nqrD (UniProtKB No. A0KHC9) was expressed.
Species	<i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.