



L-Norephedrine [HRP] (DAGA-350H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Phenylpropanolamine (PPA; Accutrim; β -hydroxyamphetamine), also known as the stereoisomers norephedrine and norpseudoephedrine, is a psychoactive drug of the phenethylamine and amphetamine chemical classes which is used as a stimulant, decongestant.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	L-Norephedrine
----------	----------------