



# Recombinant Protein Aethia pusilla (Least auklet) MT-ND6, E.coli (DAGA-3521)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	The following tags are available.
	N-terminal His-tagged and C-terminal Myc-tagged
	N-terminal His-SUMO-tagged and C-terminal Myc-tagged
	N-terminal His-B2M-tagged and C-terminal Myc-tagged
	N-terminal His-GST-tagged and C-terminal Myc-tagged
	N-terminal His-tagged
	N-terminal His-SUMO-tagged
	N-terminal His-B2M-tagged
	N-terminal His-GST-tagged
	N-terminal GST-tagged
	Tag-Free
<b>Antigen Description</b>	Recombinant 6xHis tagged Aethia pusilla (Least auklet) MT-ND6 (UniProtKB No. P43190) was expressed.
<b>Species</b>	Aethia pusilla
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer,50% glycerol
<b>Storage</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.

---