



Recombinant Protein Aethia pusilla (Least auklet) MT-ND6, Mammalian cell (DAGA-3522)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p>
Antigen Description	Recombinant 6xHis tagged Aethia pusilla (Least auklet) MT-ND6 (UniProtKB No. P43190) was expressed.
Species	Aethia pusilla
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.