



Recombinant Protein *Aethionema cordifolium* (Lebanon stonecress) atpI, E.coli (DAGA-3526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged
 N-terminal His-SUMO-tagged and C-terminal Myc-tagged
 N-terminal His-B2M-tagged and C-terminal Myc-tagged
 N-terminal His-GST-tagged and C-terminal Myc-tagged
 N-terminal His-tagged
 N-terminal His-SUMO-tagged
 N-terminal His-B2M-tagged
 N-terminal His-GST-tagged
 N-terminal GST-tagged
 Tag-Free

Antigen Description	Recombinant 6xHis tagged <i>Aethionema cordifolium</i> (Lebanon stonecress) atpI (UniProtKB No. A4QJA3) was expressed.
Species	<i>Aethionema cordifolium</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.
