



Recombinant Protein *Aethionema grandiflorum* (Persian stone-cress) psbl, E.coli (DAGA-3529)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview

The following tags are available.

- N-terminal His-tagged and C-terminal Myc-tagged
- N-terminal His-SUMO-tagged and C-terminal Myc-tagged
- N-terminal His-B2M-tagged and C-terminal Myc-tagged
- N-terminal His-GST-tagged and C-terminal Myc-tagged
- N-terminal His-tagged
- N-terminal His-SUMO-tagged
- N-terminal His-B2M-tagged
- N-terminal His-GST-tagged
- N-terminal GST-tagged
- Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

Antigen Description	Recombinant 6xHis tagged <i>Aethionema grandiflorum</i> (Persian stone-cress) psbl (UniProtKB No. A4QJI3) was expressed.
Species	<i>Aethionema grandiflorum</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg

Buffer	Tris-based buffer,50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.