



Mycophenolic acid [BSA] (DAGA-352B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Mycophenolic Acid (MPA) is an immunosuppressive agent with antitumor and antiviral effects. It inhibits the synthesis of guanosine nucleotides by binding the enzyme inosine monophosphate dehydrogenase (IMPDH). This limits the proliferation of T- and B-lymphocytes.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Mycophenolic acid; mycophenolate; MPA; MYCOPHENOLIC ACID;