



N-acetylprocainamide [KLH] (DAGA-353K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	N-Acetylprocainamide (a major metabolite of Procainamide) may cause anti-arrhythmic action, leading to cardiac toxicity and kidney failure. Procainamide is used to correct irregular heartbeats to a normal rhythm and to slow an overactive heart.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	4-(acetylamino)-n-(2-(diethylamino)ethyl)-benzamid; 4-(acetylamino)-n-[2-(diethylamino)ethyl]-benzamid; 4'-((2-(diethylamino)ethyl)carbamoyl)-acetanilid; 4'-((2-(diethylamino)ethyl)carbamoyl)acetanilide; acecainide
-----------------	--