



Recombinant Protein African swine fever virus (isolate Tick/South Africa/Pretoriuskop Pr4/1996) Pret-144, Mammalian cell (DAGA-3542)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged African swine fever virus (isolate Tick/South Africa/Pretoriuskop Pr4/1996) Pret-144 (UniProtKB No. P0CAB1) was expressed.
Species	African swine fever virus
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.