



# Recombinant Agkistrodon contortrix contortrix Zinc metalloproteinase/disintegrin, Mammalian cell (DAGA-3564)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The following tags are available.

N-terminal His-tagged and C-terminal Myc-tagged

N-terminal His-tagged

Tag-Free

The tag type will be determined during production process. If you have specified tag type, please tell us and we will develop the specified tag preferentially. If we fail in developing the specified tag, we will send other tag type directly without negotiation.

### Antigen Description

Recombinant 6xHis tagged Agkistrodon contortrix contortrix (Southern copperhead) (UniProtKB No. Q9IAB0) was expressed.

### Species

Agkistrodon contortrix contortrix

### Purity

>85% (SDS-PAGE)

### Conjugate

His

### Format

Liquid

### Concentration

Batch dependent - please inquire should you have specific requirements.

### Size

1 mg

### Buffer

Tris-based buffer, 50% glycerol

### Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.