



Recombinant Protein Agkistrodon piscivorus piscivorus (Eastern cottonmouth) MT-ND4, E.coli (DAGA-3568)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-SUMO-tagged and C-terminal Myc-tagged N-terminal His-B2M-tagged and C-terminal Myc-tagged N-terminal His-GST-tagged and C-terminal Myc-tagged N-terminal His-tagged N-terminal His-SUMO-tagged N-terminal His-B2M-tagged N-terminal His-GST-tagged N-terminal GST-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged Agkistrodon piscivorus piscivorus (Eastern cottonmouth) MT-ND4 (UniProtKB No. P92503) was expressed.
Species	Agkistrodon piscivorus piscivorus
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol

Storage

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
