



Phenytoin [HRP] (DAGA-356H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Phenytoin is an anticonvulsant drug which can be useful in the treatment of epilepsy. Phenytoin blocks voltage-sensitive sodium channels in neurons. This action leads to a delay in neuronal electrical recovery. |
| Conjugate | HRP |
| Applications | ELISA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

Keywords Dilantin; Diphenylhydantoin; Epanutin; Phenytek;