



Phenytoin [KLH] (DAGA-356K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Phenytoin is an anticonvulsant drug which can be useful in the treatment of epilepsy. Phenytoin blocks voltage-sensitive sodium channels in neurons. This action leads to a delay in neuronal electrical recovery.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Dilantin; Diphenylhydantoin; Epanutin; Phenytek;
-----------------	--