



# Recombinant Protein Agrobacterium tumefaciens (strain C58) kup2, Mammalian cell (DAGA-3577)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | <p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p> |
| <b>Antigen Description</b> | Recombinant 6xHis tagged Agrobacterium tumefaciens (strain C58) kup2 (UniProtKB No. Q8UJL0) was expressed.                                   |
| <b>Species</b>             | Agrobacterium tumefaciens                                                                                                                    |
| <b>Purity</b>              | >85% (SDS-PAGE)                                                                                                                              |
| <b>Conjugate</b>           | His                                                                                                                                          |
| <b>Format</b>              | Liquid                                                                                                                                       |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.                                                                      |
| <b>Size</b>                | 1 mg                                                                                                                                         |
| <b>Buffer</b>              | Tris-based buffer, 50% glycerol                                                                                                              |
| <b>Storage</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                          |