



# Recombinant Protein Agrobacterium radiobacter (strain K84 / ATCC BAA-868) cobS, Yeast (DAGA-3578)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	The following tags are available. N-terminal His-tagged Tag-Free
<b>Antigen Description</b>	Recombinant 6xHis tagged Agrobacterium radiobacter (strain K84 / ATCC BAA-868) cobS (UniProtKB No. B9JG16) was expressed.
<b>Species</b>	Agrobacterium radiobacter
<b>Purity</b>	>85% (SDS-PAGE)
<b>Conjugate</b>	His
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.