



# Recombinant Protein Agrobacterium radiobacter (strain K84 / ATCC BAA-868) atpB, Baculovirus (DAGA-3580)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |  |
|----------------------------|--|
| <b>Product Overview</b>    | <p>The following tags are available.</p> <p>N-terminal His-tagged and C-terminal Myc-tagged</p> <p>N-terminal His-tagged</p> <p>Tag-Free</p> |
| <b>Antigen Description</b> | Recombinant 6xHis tagged Agrobacterium radiobacter (strain K84 / ATCC BAA-868) atpB (UniProtKB No. B9JA28) was expressed.                    |
| <b>Species</b>             | Agrobacterium radiobacter  |
| <b>Purity</b>              | >85% (SDS-PAGE)  |
| <b>Conjugate</b>           | His  |
| <b>Format</b>              | Liquid   |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.  |
| <b>Size</b>                | 1 mg   |
| <b>Buffer</b>              | Tris-based buffer, 50% glycerol  |
| <b>Storage</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |