



Recombinant Protein Agrobacterium rhizogenes kup, E.coli (DAGA-3583)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-SUMO-tagged and C-terminal Myc-tagged N-terminal His-B2M-tagged and C-terminal Myc-tagged N-terminal His-GST-tagged and C-terminal Myc-tagged N-terminal His-tagged N-terminal His-SUMO-tagged N-terminal His-B2M-tagged N-terminal His-GST-tagged N-terminal GST-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged Agrobacterium rhizogenes kup (UniProtKB No. Q9F5A5) was expressed.
Species	Agrobacterium rhizogenes
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer,50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to

one week.
